powertube®

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
Instruction manual

This instruction manual is considered an integral part of the device. It shall be kept within reach at all times. Exact compliance with these instructions is essential for the correct handling of the device as well as for the safety of patient and user.

Tabel of contents Page

1.	Illustration of the device	4
2.	Control panel with functions	5
3.	Description of the device	6
4.	Quick start guide	6
5.	Battery and battery compartment	7
6.	Specifications	8
7.	Handling	9
8.	Indications	10
9.	Contraindications	10
10.	Safety tips	10
11.	Care, maintenance and storage	11
12.	Errors and malfunctions, function test	12
13.	Accessories	12
14.	Guarantee	12
15.	Safety inspection	13
16.	Conformity declaration CE 0482	14
I.	Risk of mutual interference see annex «Technical description EN/IEC 60601-1-2»	

1. Illustration of the device

Therapy head

Plug connection for cable to handhold

Insulation

Handhold with cable can be used instead of therapy head


Control panel and display

Gold- / Silver-plated casing with integrated electronics

Cable with plug for connection to PowerTube® Gold/Silver

Cover of battery compartment. Use a coin to open. Requires one IEC 6LR61 9-volt alkaline battery.

2. Control panel with functions

-  1. Turning on the device (push button briefly)
-  2. Select number of minutes (1 - 7 min.)
-  3. Start therapy (push button briefly) (manual STEP selection)
-  4. Turning off the device (push button for an extended period of time or device will automatically turn off after 30 seconds)

3. Description of the device

Intended use

The PowerTube® Gold/Silver is intended for the electronic stimulation of nerves. It is used for self-treatment of symptoms and pain.

Its application does not replace treatment or instructions from a doctor, but has a positive influence on symptoms and pain.

Using the device can alleviate pain. Please refer to the instruction manual for detailed information on the «intended use».

Function designation

The transcutaneous electronic nerve stimulation device PowerTube® Gold/Silver is an easy to use, battery-operated electronic device designed to stimulate the nerves.

Corresponding nerve points are stimulated directly with electronic impulses. The device can be operated by persons who do not have formal medical training.

The device is described in the enclosed easy to follow instruction manual to make self-treatment a straightforward process. The handling of the device is intended to suit users with no medical training. Various applications as well as corresponding therapy times are outlined in the instruction manual.

4. Quick start guide

Before using for the first time

Immediately after receipt, the device should be inspected for transport damages and be subject to a function test. In case of transport damage, the shipper shall be notified immediately and a damage report be drawn up to support any claim for damages. In addition to name and address of the recipient, the protocol shall include the date of receipt, type and serial number of the device as well as a description of the damage. The original packaging should be kept during the guarantee period to ensure that the device can be returned in its original packaging if necessary.


Every device is tested by the manufacturer in terms of functionality and safety before being shipped. Please observe the following to ensure that the device continues to operate properly following transport and installation:

The operator should only start using the device after carrying out a function test (see chapter 12).

5. Battery and battery compartment

Use alkaline batteries only and check the polarities (plus / minus) carefully. By using this kind of battery you may be certain to end your therapy time even in case the battery is almost empty by showing the sign «low battery». The device is adjusted to this type of batteries.

If you use rechargeable accus (e. g. NiMH), the proper functioning of the sign «low battery» is not completely reliable because of the different mode of discharging. The device may abruptly stop functioning without prior warning. If you still prefer to use rechargeable accus, please recharge them periodically.

6. Specifications

Display	<ul style="list-style-type: none"> • Time of therapy in minutes • 7 light-emitting diodes • Therapy steps 1 - 3 • Exchange light-emitting diode battery
Error messages	<p>Intermittent tone: handholds or contact areas are not applied properly</p> <p>Continuous tone: handholds have short circuited</p>
Battery	IEC 6LR61 9-volt block battery
Output voltage	AC max. 9-volt pp / min. 6-volt pp
Duration of treatment	3 - 21 minutes
Programme sequence	3 steps
Measurements	Length 18 cm, diameter 3,1 cm
Weight	~ 1100 g
Application class	BF
Safety class	Ila unit
IP protection class	<p>Device: IP40 Protection against entry of solid foreign substances $\varnothing > 1$ mm, not protected against water.</p> <p>Softbag / Case: IPX2 Protection against dripping water (tilted) up to 15° relative to vertical line.</p>
Operation	<ul style="list-style-type: none"> • 5 °C - 40 °C • humidity 15 % - 93 %, without cond. • atm. pressure 700 hPa - 1060 hPa
Transport and storage	<ul style="list-style-type: none"> • -25 °C - 70 °C • humidity - 93 %, without cond. • atm. pressure 700 hPa - 1060 hPa
Graphic symbols	 Please refer to instruction manual!
	 Typ BF
Production date	month/year-serial number

7. Handling

1. Turning on the device: push I / GO button briefly.
2. Set appropriate therapy time per step by means of TIME button. Start with 1 minute per step, i. e. the complete treatment time is 3 minutes. Increase the time with every treatment by 1 minute per step. Complete steps 1 and 2 (see below and beside) once a day.
3. Place device on trigger point.
4. Start treatment by pushing the I / GO button.
5. You will hear a beep signal after completion of every step.
6. The end of the treatment is announced with 2 beep signals.
7. Turning off the device: push I / GO button for a prolonged periode of time or device turns off automatically after 30 seconds.

Error message:

Insufficient or no skin contact is indicated with continuous beep signals.

Step 1

Local treatment at trigger point

Step 2

Treatment of abdominal area

Please ensure good body contact. Moisten skin with water if necessary.

Consume large quantities of water before and after the treatment.

8. Indications

- Pain relief

9. Contraindications

- Pregnancy, cardiac pacemaker, open wounds and fresh scars.
- Do not use in combination with other devices.

10. Safety tips

- While operating in immediate vicinity of a short- or microwave may lead to instability of the stimulator-output.
- PowerTube® Gold/Silver may only be used for the electrical stimulation of nerves and musculature.
- Do not use the device for patients who have been fitted with a cardiac pacemaker.
- Please note that risk of strangulation by cable because of the excess length.
- Therapie times especially first use of the device (without cable, direct therapy):
 - Head max. 2-2-2 minutes
 - Thorax max. 3-3-3 minutes
- After one week, increase slightly the time of therapy.

- Simultaneous use with other devices is prohibited.
- PowerTube® Gold/Silver may only be operated with a 9-volt battery.
- PowerTube® Gold/Silver may only be operated with the accessories included in the range of delivery.
- Famedtec GmbH is only responsible for the reliability and safety of PowerTube® Gold/Silver if:
 - PowerTube® Gold/Silver is operated in accordance with the instruction manual,
 - expansions and modifications are carried out only by persons who are authorised by Famedtec GmbH.
- Please recycle the batteries in accordance with the national regulations.
- The expected life of the device, including parts and accessories amounts 10 years.
- PowerTube® Gold/Silver should be disposed of as hazardous waste.
- For technical questions in connection with the handling, or maintenance contact the manufacturer.

11.Care, maintenance and storage

Effective protection of the PowerTube® Gold/Silver from damage includes:

- Appropriate operation and maintenance.
- Protection from moisture, high temperatures, dust and direct sunlight.
- Protection from children, animals and pests.
- Wipe electrodes after every use to disinfect them.
- Remove the battery if the device is not being used for a prolonged period of time.
- The expected life of new alkaline batteries amounts 2 months (daily duration of 20 minutes).

12. Errors and malfunctions, function test

- Discharge of the battery is displayed on a battery indicator LED.
- After being turned on, PowerTube® Gold/Silver automatically runs a battery test.
- An acoustic warning signal during therapy indicates that the electrodes have not been properly applied.

13. Accessories

- Hand electrode
- Instruction manual
- IEC 6LR61 9-volt battery
- Softbag / Case for PowerTube® Gold/Silver

14. Guarantee

Device and accessories shall be inspected for transport damage immediately after receipt. Corresponding claims for damages can only be asserted if the vendor or shipper have been notified immediately. A damage report shall be compiled.

- The guarantee for PowerTube® Gold/Silver is two years from day of delivery.
- The extent of guarantee includes free repair of the device, provided the damage was caused by a material or manufacturing error.

- Additional claims, particularly claims for damages, are excluded.
- The repair may only be carried out by Famedtec GmbH, one of its agencies or authorised specialist dealers. The guarantee shall become void if improper modifications or repairs are carried out.
- The guarantee period will neither be extended nor renewed because of guarantee work.
- Famedtec GmbH reserves all rights to this instruction manual, particularly the right of duplication and distribution as well translation. No part of this instruction manual may be copied or distributed in any form without prior written consent from Famedtec GmbH.
- Information contained in this instruction manual may be modified or upgraded by Famedtec GmbH (updates can be found on www.famedtec.ch).

15. Safety inspection

§ 6 of the directive on installing, operating and using active medical products stipulates regular safety inspections to guarantee the safety of the device against age, wear and tear etc.

The operator is responsible for having these safety inspections carried out on the device to the prescribed extent and in due time. Safety inspections may only be carried out by the manufacturer or by explicitly authorised persons.

The following safety inspections are required for Power-Tube® Gold/Silver:

- Examination of labelling and instruction manual
- Visual inspection of device and accessories for damages
- Function test of all sensors and control lamps on the device

16. Conformity declaration

CE 0482

Name and address of the company	Famedtec GmbH Hauptstrasse 42 CH-3255 Rapperswil BE Switzerland
We declare in the position as solely responsible company that:	
the medical electrical equipments	PowerTube® Gold Typ FPT.B1 PowerTube® Silver Typ FPT-D.B2 according to appendix IX of the directives 93/42/EWG
in safety class	Ila according to appendix IX of the directives 93/42/CEE
all required standards of the directives mentioned in the medical dispositions 93/42/EWG are herewith completely fulfilled and adhered to.	
Conformity evaluation procedure	Directives 93/42/EWG appendix V
Registration no.	7288DE414150310
Conformity of the quality management system is acknowledged and monitored by	MEDCERT GmbH Pilatuspool 2 D-20355 Hamburg Deutschland CE 0482
CH-Rapperswil BE, 2017	Fritz Wermuth (CEO Famedtec GmbH)
The period of validity of the certificate	23 January 2020

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Instruction manual
PowerTube®

© 2017 Famedtec GmbH. All rights reserved. All trade names are registered trademarks of

Famedtec GmbH
Hauptstrasse 42
CH-3255 Rapperswil BE
Switzerland

Edition 2017-04-01
Printed in Switzerland

powertube®
Powertube TENS Therapie Technik